MDS Lite User Guide

CE

MDS Lite User Guide

MDS-Lite-Benutzerhandbuch

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Guide d'utilisation de MDS Lite

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English

Deutsch

Español

Français

Italiano

Português

SonoSite, Inc.

21919 30th Drive SE

Bothell, WA 98021

USA

T: 1-888-482-9449 or 1-425-951-1200

F: 1-425-951-1201

SonoSite Ltd

Alexander House

40A Wilbury Way

Hitchin

Herts SG4 OAP

UK

T: +44-1462-444800

F: +44-1462-444801

Caution: | United States federal law restricts this device to sale by or on the order of a physician.

"MicroMaxx," "TITAN," and "SonoSite TITAN" are trademarks or registered trademarks of SonoSite, Inc.

Non-SonoSite product names may be trademarks or registered trademarks of their respective owners.

SonoSite products may be covered by one or more of the following U.S. patents: 4454884, 4462408, 4469106, 4474184, 4475376, 4515017, 4534357, 4542653, 4543960, 4552607, 4561807, 4566035, 4567895, 4581636, 4591355, 4603702, 4607642, 4644795, 4670339, 4773140, 4817618, 4883059, 4887306, 5016641, 5050610, 5095910, 5099847, 5123415, 5158088, 5197477, 5207225, 5215094, 5226420, 5226422, 5233994, 5255682, 5275167, 5287753, 5305756, 5353354, 5365929, 5381795, 5386830, 5390674, 5402793, 5,423,220, 5438994, 5450851, 5456257, 5471989, 5471990, 5474073, 5476097, 5479930, 5482045, 5482047, 5485842, 5492134, 5517994, 5529070, 5546946, 5555887, 5603323, 5606972, 5617863, 5634465, 5634466, 5636631, 5645066, 5648942, 5669385, 5706819, 5715823, 5718229, 5720291, 5722412, 5752517, 5762067, 5782769, 5800356, 5817024, 5833613, 5846200, 5860924, 5893363, 5916168, 5951478, 6036643, 6102863, 6104126, 6113547, 6117085, 6142946, 6203498 B1, 6371918, 6135961, 6364839, 6383139, 6416475, 6447451, 6471651, 6569101, 6575908, 6604630, 6648826, 6835177, 6817982, 6730035, 6962566, 7169108, D0280762, D0285484, D0286325, D0300241, D0306343, D0328095, D0369307, D0379231, D456509, D461895, D509900. Other patents pending.

P05401-03 07/2007

Copyright 2007 by SonoSite, Inc.

All rights reserved.

MDS Lite User Guide

Contents

Introduction	1
Safety	2
MDS Lite Features	4
Connecting Ultrasound System	8
Connecting Printer and DVD Recorder or VCR	8
Connecting Triple Transducer Connect	14
Cables, Connectors, and Mini-Dock Symbols	15
Troubleshooting	19
Cleaning and Disinfecting the MDS Lite	20
Cleaning the Peripherals	20
Specifications	21

Introduction

The Mobile Docking System Lite (MDS Lite) provides a mobile work platform for transducers, accessories, peripherals, and one of the following ultrasound systems:

- MicroMaxx® ultrasound system
- TITAN® high-resolution ultrasound system

This user guide provides information on the MDS Lite and on connecting accessories and peripherals. For information on the peripherals, please see the manufacturer's operating instructions shipped with each product. For information on the ultrasound system and transducers, see the ultrasound system user guide.

WARNING: To avoid the risk of patient injury, misdiagnosis, or user injury, read all the warnings in the ultrasound system user guide and supplemental guides.

Conventions

These conventions are used in this user guide:

- A **WARNING** describes precautions necessary to prevent injury or loss of life.
- A **Caution** describes precautions necessary to protect the products.
- Numbered steps in procedures must be performed in a specific order.
- Bulleted lists present information in list format but do not imply a sequence.

Labeling Symbols and Terms

Labeling symbols and terms used on the system are explained in the ultrasound system user guide.

Customer Assistance and Comments

Questions and comments are encouraged. SonoSite is interested in your feedback regarding the system and the user guide. Please call SonoSite at **1-888-482-9449**. For technical support, call **1-877-657-8118** or if you are outside the USA, call the nearest SonoSite representative. You can also e-mail SonoSite at **comments@sonosite.com**.

Safety

For maximum safety, observe the following warnings and cautions before using the accessories and peripherals with the ultrasound system and transducers.

Electrical Safety

- WARNING:** To avoid the risk of electrical shock and fire hazard, a multiple portable socket outlet (MPSO) or extension cord should not be connected to the auxiliary mains socket outlets on the MDS Lite.
- To avoid the risk of electrical shock, use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade" or the equivalent. The grounding wire must not be removed or defeated.
- To avoid the risk of electrical shock and fire hazard, inspect the power supply, plugs, cables, and power cords on a regular basis for damage.
- To avoid the risk of electrical shock, use only accessories and peripherals recommended by SonoSite, including the power supply. Connection of accessories and peripherals not recommended by SonoSite could result in electrical shock. Contact SonoSite or your local representative for a list of accessories and peripherals available from or recommend by SonoSite.
- To avoid the risk of electrical shock, use commercial grade peripherals recommended by SonoSite on battery power only. Do not connect these products to AC mains power when using the system to scan or diagnose a patient/subject. Contact SonoSite or your local representative for a list of the commercial grade peripherals available from or recommended by SonoSite.
- To avoid possible electrical shock or electromagnetic interference, verify proper operation and compliance with relevant safety standards for all equipment before clinical use. Connecting additional equipment to the ultrasound system constitutes configuring a medical system. SonoSite recommends verifying that the system, all combinations of equipment, and accessories connected to the ultrasound system comply with JACHO installation requirements and/or safety standards such as AAMI-ES1, NFPA 99 OR IEC Standard 60601-1-1 and electromagnetic compatibility standard IEC 60601-1-2 (Electromagnetic compatibility), and are certified according to IEC Standard 60950 (Information Technology Equipment (ITE)).

Equipment Safety

- Caution:** Excessive bending or twisting of cables can cause a failure or intermittent operation.
- Improper cleaning or disinfecting of any part of the accessories and peripherals can cause permanent damage. For cleaning and disinfecting instructions, see [“Cleaning and Disinfecting the MDS Lite”](#) on page 20 and the manufacturer’s operating instructions.
- Accessible metal of the mini-dock is not protectively earthed. Do not perform high current ground impedance test involving this part.

Biological Safety

- WARNING:** Non-medical (commercial) grade peripheral monitors have not been verified or validated by SonoSite as being suitable for diagnosis.

Electromagnetic Compatibility (EMC)

See the safety section of the ultrasound system user guide and supplemental guides for information on electromagnetic compatibility with IEC 60601-1-2, including the Manufacturer’s Declaration.

MDS Lite Features

The MDS Lite provides a mobile work platform and provides a storage area for the transducers and other supplies. The MDS Lite includes the mini-dock with connections to accessories and peripherals.

WARNING: To avoid the risk of electrical shock or injury, the equipment shall be connected to a center-tapped single phase supply circuit when users in the United States connect the equipment to a 240V supply system.


Figure 1 MDS Lite: Front and Back Views


Table 1: MDS Lite Features

Number	Description
1	System platform to hold transducer, connector, gel, and other items
2	Handles to transport and hooks for transducer cables
3	Shelf for black-and-white printer or storage
4	Area for VCR, DVD recorder, color printer, or storage

Table 1: MDS Lite Features (Continued)

Number	Description
5	Locking wheels
6	System latches
7	Security straps
8	Column
9	Height adjustment knob for system platform
10	Power cord hooks
11	Power strip
12	Security strap slots

- Raise or lower shelf**
- 1 Loosen the two bolts under the shelf with a T-handle hex wrench.
 - 2 Raise or lower the shelf to desired position.
 - 3 Tighten the two bolts under the shelf using a T-handle hex wrench until there is no gap between the shelf and column. See [Figure 2](#).
 - 4 Grasp the sides of the shelf and lift up and then push down to ensure that it is securely attached.

**Figure 2 Tighten Shelf Bolts**

-
- | | | |
|---------------------------------------|---|---|
| Raise or lower system platform | 1 | Loosen the adjustment knob. See Figure 3 . |
| | 2 | Raise or lower the system platform to desired position. |
| | 3 | Tighten the adjustment knob until system platform is flush against the column and securely tightened. |
| | 4 | Grasp the sides of the system platform and lift up and then push down to ensure that it is securely attached. |


Figure 3 Adjustment Knob

-
- | | | |
|----------------------|---|--|
| Insert system | 1 | Rotate system latches outward. |
| | 2 | Place the ultrasound system on the system platform, ensuring that the cable is positioned between the handles of the MDS Lite. |
| | 3 | Secure by rotating system latches inward. |

Caution:	To avoid damaging the transducer cable, place transducer cables between front handles.
-----------------	--

-
- | | | |
|----------------------|---|-------------------------------|
| Remove system | 1 | Rotate latches outward. |
| | 2 | Remove the ultrasound system. |

Attach mini-dock

- 1 Turn off the ultrasound system.
The mini-dock may not function properly if you connect the system while the power is on.
- 2 Place the mini-dock on the system platform.
- 3 Slide the mini-dock feet under the ultrasound system and firmly press the mini-dock into the system connector.

The mini-dock is correctly attached when there is no gap between the system and mini-dock. See [Figure 4](#).

Note: If at any time the ultrasound system should lose connectivity, turn off the system and reattach the mini-dock and system.

Note: Always verify that the mini-dock and system are correctly attached after changing a transducer. When you change a transducer, the mini-dock may loosen.


Figure 4 Attach Mini-Dock

Lock or unlock wheels Press down on the locking levers on the front wheels to lock, and push up to unlock.

Connecting Ultrasound System


- | | |
|-------------------------|---|
| Connect AC power | <ol style="list-style-type: none">1 Ensure that the mini-dock is connected to the system.2 Connect the cables. See Figure 5.<ul style="list-style-type: none">• DC power cord (A)• AC power cord (I)3 If connecting system only, connect the system AC power cord (K) to a hospital-grade electrical outlet. |
|-------------------------|---|

Connecting Printer and DVD Recorder or VCR

WARNING:	To avoid the risk of electrical shock or injury, connect only accessories or peripherals to the MDS Lite that meet the following requirements: • Does not exceed two accessories or peripherals • Is not the color and black-and-white printers
-----------------	---

Caution:	To avoid damaging the system, use only peripherals recommended by SonoSite with the ultrasound system.
-----------------	--

The MDS Lite and mini-dock provide connections for a printer and a DVD recorder or VCR. See [Figure 5](#). For specific information on using the DVD recorders, VCR, and printers, including warnings and cautions, see the manufacturer's operating instructions.

- * If Sony VCR, set dip switches: 1-4 Down and 5,6 Up.
If Panasonic VCR, set switch to Line.
- ** If Sony DVD recorder, use adapter G2 with E.
If Panasonic DVD recorder, use adapter G with E.
If Panasonic VCR, use adapter F with E.

Figure 5 Connectivity of System, Printer, and a DVD Recorder or VCR

DVD Recorder

Connect cables	1	Lock the wheels.
	2	Insert security strap in middle slots, if necessary.
	3	Move the device forward to allow access to the connectors, if necessary.
	4	Ensure that the mini-dock is connected to the system.
	5	Connect the cables. See Figure 5 on page 9 and “ Cables, Connectors, and Mini-Dock Symbols ” on page 15.
		• S-video cables (B) (B). Insert the ferrite end of the cable into the mini-dock. • RS-232 cable (E): • For Sony, attach DVD adapter (G2) • For Panasonic, attach DVD adapter cable (G) • AC power cord (I)
	6	Connect system power cord (K) to a hospital-grade electrical outlet.
	7	Fold excess cable behind the device, and insert the cables on either side of column, using the grommets to secure cables, if necessary.
	8	Tighten the strap.

Set up system	1	Press the Setup key on the ultrasound system.
	2	Select Connectivity .
	3	In the Video Mode list, select the appropriate video format: NTSC or PAL .
	4	In the Serial Port list, select DVD .
	5	Restart the system to activate DVD connectivity.

Note: Because the bar code scanner, DVD recorder, and VCR require the same port on the mini-dock, you can connect only one of these peripherals at a time.

Set up DVD recorder for recording (Panasonic)	1	Turn on the DVD recorder.
	2	Press the Function key (f) then V on the system.
	3	Press Shift and Menu on the DVD recorder. The DVD menu appears on the system.
	4	Select Setup and then press Enter . A help diagram is displayed on the screen indicating how to use the DVD buttons.
	5	From the Remote menu, select Baud Rate and 9600 .
	6	From the Remote menu, select Protocol and VTR1 .
	7	From the Picture menu, select Input Select and S-Video .
	8	From the Connection menu, select TV System and NTSC or PAL .
	9	Press Return on the DVD recorder to exit.
	10	On the system, press the Function key (f) then V to exit.

-
- Record images**
- 1 Turn on the DVD recorder.
 - 2 Insert a DVD.
 - 3 Press the **Record** key on the ultrasound system, or use controls on the DVD recorder to start and stop recording.
Each time the Record key is pressed, a new file is saved to the DVD.
(Panasonic Only) The DVD can be viewed on a PC if finalized on the DVD recorder.
-

- Play back images**
- Note:* Do a test recording before a patient exam to ensure that all connectivity is correct.
- 1 Press the Function key (**f**) then **V** on the system.
 - 2 Press **Play** on the DVD recorder. Follow the manufacturer's instructions for playing back images.

Connect cables	1 Lock the wheels. 2 Insert security strap in outer slots, if necessary. 3 Move the device forward to allow access to the connectors, if necessary. 4 Ensure that the mini-dock is connected to the system. 5 Connect the cables. See Figure 5 on page 9 and “Cables, Connectors, and Mini-Dock Symbols” on page 15. • S-video cables (B) (B). Insert the ferrite end of the cable into the mini-dock. • RS-232 cable (E), including audio left/right input/output. For Panasonic VCR, attach VCR adapter (F) to E. • AC power cord (I) 6 Connect system power cord (K) to a hospital-grade electrical outlet. 7 Fold excess cable behind the device, and insert the cables on either side of column, using the grommets to secure cables, if necessary. 8 Tighten the strap.
Set up system	1 Press the Setup key on the ultrasound system. 2 Select Connectivity. 3 In the Video Mode list, select the appropriate video format: NTSC or PAL. 4 In the Serial Port list, select VCR. 5 Restart the system to activate VCR connectivity. <i>Note:</i> Because the bar code scanner, DVD recorder, and VCR require the same port on the mini-dock, you can connect only one of these peripherals at a time.
Record images	1 Ensure that the cables are correctly installed. 2 Turn on the VCR. 3 Ensure that the input switch is correctly set. (See the manufacturer’s instructions for location of switches.) • Sony: Dip switches: 1-4 Down and 5, 6 Up • Panasonic: Line 4 Press the Record key or use the controls on the VCR to start and stop recording.
Play back images	<i>Note:</i> Do a test recording before a patient exam to ensure that all connectivity is correct. 1 Turn on the VCR. 2 Press the Function key (f) then V on the system. 3 Press Play on VCR.

Printer

Connect cables	1	Lock the wheels.
	2	Move the printer forward to allow access to the connectors, if necessary.
	3	Ensure that the mini-dock is connected to the system.
	4	Connect the cables. See Figure 5 and “Cables, Connectors, and Mini-Dock Symbols” on page 15. • Print control cable (H). Insert the ferrite end of the cable into the mini-dock. • Printer AC power cord (J) • Composite video cable (C)
	5	Connect system power cord (K) to a hospital-grade electrical outlet.
	6	Insert the cables in either side of the column, folding excess cable behind the device, and using the grommets to secure cables, if necessary.
	7	Tighten the strap.
Set up system	1	Press the Setup key on the ultrasound system.
	2	Select Connectivity .
	3	In the Printer list, select the desired printer.
		Only the printers listed are recommended for use with the system.
Print images		<i>Note:</i> Do a test printing before a patient exam to ensure that all connectivity is correct.
	1	Ensure that the cables are correctly installed.
	2	Turn on the printer.
	3	(Color Printer Only) Change input setting from S-video to video.
	4	Press the Print key or use the controls directly on the printer.

Connecting Triple Transducer Connect

The Triple Transducer Connect is a three-transducer module for the MDS Lite that allows connection of up to three transducers simultaneously while the ultrasound system is docked.

Connect transducer connector to system	1	Rotate system latches outward.
	2	Remove any transducer attached to the system.
	3	Pull the transducer latch up and rotate 90°.
	4	Hold the system at an angle.
	5	Align and insert the transducer connector to the connector on the bottom of the ultrasound system.
	6	Turn the latch back 90° and press in, securing the transducer connector to the ultrasound system.
	7	Place the ultrasound system on the system platform, ensuring that the cable is positioned between the handles of the MDS Lite.
	8	Secure by rotating system latches inward.


Figure 6 Connect Transducer Connector to System

Connect transducer to Triple Transducer Connect	1	Pull the transducer latch up and rotate 90°.
	2	Align and insert the transducer connector to the connector on the bottom of the Triple Transducer Connect.
	3	Turn the latch and press down, securing the transducer connector to the Triple Transducer Connect.

Select transducer

Note: When you use the Triple Transducer Connect, there may be a minor reduction of transducer penetration.

Push the transducer select button on the Triple Transducer Connect to select the desired transducer. See [Figure 6](#) on page 14.

- When selecting another transducer, the system restarts. The light on the button flashes momentarily and then turns to a solid green.
- Only one transducer is accessible at a time.

Cables, Connectors, and Mini-Dock Symbols

The cable kit includes cables for both the MDSe and the MDS Lite. You need to use only the cables for your product.

Table 2: Connectivity Symbols on Mini-Dock

Symbol	Definition
	DC input
	Print control
	USB
	Ethernet
	RS-232 (DVD recorder, VCR, or bar code scanner)
	S-video out
	S-video in
	DVI video out
	Composite video out
	Audio out
ECG	ECG/Footswitch

Table 3: Manually Installed Cables and Connectors










Letter	Name	Illustration	Connector 1 [[Illustration] [Connects to]	Connector 2 [[Illustration] [Connects to]	Product Used On
B	S-video cable (playback) 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
B	S-video cable (record) 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
B1	S-video extension cable (Sony DVD recorder only) 12 in./30.5 cm		 S-video adapter	 DVD recorder	MDSe
B2	S-video adapter (Sony DVD recorder only)		 S-video cable	 S-video extension cable	MDSe
C	Composite video cable 6 ft./1.8 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
D	Composite video cable 12 in./30.5 cm		 Printer	 DVD recorder or VCR	MDSe

Table 3: Manually Installed Cables and Connectors





















Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
E	RS-232 cable (control/audio) 6.5 ft./2 m		 Mini-dock	 DVD recorder or VCR	MDS Lite
E1	Audio extension cables (Sony DVD recorder only) 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSe
E1	Audio extension cables (Sony DVD recorder only) 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSe
F	VCR adapter (Panasonic only) 12 in./30.5 cm		 RS-232 cable	 VCR	MDSe MDS Lite
G	DVD adapter (Panasonic only) 12 in./30.5 cm		 RS-232 cable	 DVD recorder	MDSe MDS Lite

Table 3: Manually Installed Cables and Connectors

Letter	Name	Illustration	Connector 1 [[Illustration] [Connects to]	Connector 2 [[Illustration] [Connects to]	Product Used On
G2	DVD adapter (Sony only) 12 in./30.5 cm		 RS-232 cable	 DVD recorder	MDSe MDS Lite
H	Print control cable 6 ft./1.8 m		 Mini-dock	 Printer	MDS Lite
I	AC power cord 18 in./45.7 cm		 Power supply	 Power strip	MDSe MDS Lite
I	AC power cord 18 in./45.7 cm		 DVD recorder or VCR	 Power strip	MDSe MDS Lite
I	AC power cord 18 in./45.7 cm		 Printer	 Power strip	MDSe MDS Lite
J	AC power cord 39 in./1 m		 Printer	 Power strip	MDS Lite
K	System AC power cord 10 ft./3 m		 Power strip	 Wall outlet	MDSe MDS Lite (ships in separate box)

Troubleshooting

Symptom	Solution
DVD does not record.	- Check the serial port setting in system setup, including Record key setting and DVD set up. See "Set up system" on page 10 and "Set up DVD recorder for recording (Panasonic)" on page 10. - Check the cable connections on stand and DVD. See "Connecting Printer and DVD Recorder or VCR" on page 8. - Check the DVD to ensure that it is turned on and set up properly. See the DVD manufacturer's instructions, if necessary.
DVD does not play back.	- Check the cable connections on stand and DVD. (See "Connecting Printer and DVD Recorder or VCR" on page 8. - Check the display setting to ensure that it is set to S-video. Press the +/- key to toggle video format on the display (default setting is Composite).
VCR does not record.	- Check the video format setting in system setup, including Print/VCR key setting. See "Set up system" on page 12. - Check the cable connections on stand and VCR. See "Connecting Printer and DVD Recorder or VCR" on page 8. - Check the VCR to ensure that it is turned on and set up properly. See the VCR manufacturers' instructions, if necessary. - Check the input switch is properly set. See "Set up system" on page 13.
VCR does not play back.	- Check the cable connections on stand and VCR. See "Connecting Printer and DVD Recorder or VCR" on page 8. - Check the display setting to ensure that it is set to S-Video. Press the +/- key to toggle video format on the display (default setting is Composite). - Check the switch input is properly set. See "Record images" on page 12.
Print does not work.	- Check the printer selection in system setup. See "Set up system" on page 13. - Check the cable connections. See "Connecting Printer and DVD Recorder or VCR" on page 8 or "Printer" on page 13. - Check the printer to ensure that it is turned on and set up properly. See the printer manufacturer's instructions, if necessary.

Cleaning and Disinfecting the MDS Lite

WARNING:	To avoid electrical shock, remove the system and extra battery from the mini-dock and unplug from power supply before cleaning. To avoid injury, always use protective eyewear and gloves when performing cleaning and disinfecting procedures. To avoid infection, if a pre-mixed disinfectant solution is used, observe the solution expiration date, and ensure that the date has not passed. To avoid infection, the level of disinfection required for a product is dictated by the type of tissue it contacts during use. Ensure the solution strength and duration of contact are appropriate for the equipment. For information, see the disinfectant label instructions and the recommendations of the Association for Professionals in Infection Control and Epidemiology (APIC) and FDA.
Caution:	Use only recommended cleaners or disinfectants on surfaces. Immersion-type disinfectants are not tested for use on surfaces.

Clean and disinfect	The exterior surface of the MDS Lite can be cleaned and disinfected using a recommended cleaner or disinfectant. 1 Remove the system, unplug the power supply, and detach any cables. 2 Clean the surfaces using a soft cloth lightly dampened in a mild soap or detergent cleaning solution. Apply the solution to the cloth rather than the surface. 3 Wipe the surfaces with the disinfecting solution. Theracide disinfectant is recommended. 4 Air dry or towel dry with a clean cloth.
----------------------------	---

Cleaning the Peripherals

For information about cleaning the peripherals, see the manufacturer's operating instructions.

Specifications

MDS Lite Dimensions

Width: 20 in. (50.8 cm)
Length: 30 in. (76.2 cm)
Height: 45 in. (114.3 cm)
Height range: 34-42 in. (86-107 cm)
Weight: 44 lbs. (20 kg)

Temperature and Humidity Limits

Operating: 10–40°C; 15–95% R.H.
Shipping/Storage: -35–65°C; 15–95% R.H.

Peripherals

See the manufacturer's specifications for the peripherals.

Medical Grade

Black-and-white printer

Recommended sources for printer paper: Contact Sony at **1-800-686-7669** or **www.sony.com/professional** to order supplies or to obtain the name of the local distributor.

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